

LUMECCA *peak*



PROCEDURE: Lumecca Peak

NUMBER OF TREATMENT(S): Single treatment

TIME ELAPSED: 4 weeks between Before & After image

PARAMETERS: Wavelength 515nm, Fluence 11, Cooling: Strong, Pulse Mode: Short & Long, Pulse: 789

This case highlights a female patient treated with Lumecca technology for pigmentation and overall skin clarity. Utilizing a single treatment session with strong cooling and customized pulse settings, the patient experienced a noticeable reduction in visible skin imperfections just 4 weeks post-treatment. Results included a more even skin tone and radiant complexion.

Photo Courtesy: Josie Cervantes FNP



The treatment technique is presented to demonstrate the method utilized by the practitioner. The information on the products and procedures contained in this Case Study does not represent and does not constitute medical advice or recommendations by InMode and should not be relied upon as such. It is the clinicians' responsibility to decide which application and settings they use in each clinical situation.